


AFBR-S20M2XX Qmini Sensitivity Performance



Qmini UV	Qmini VIS	Qmini NIR	Qmini WIDE-UV Sensitivity optimized for UV at 300 nm	Qmini WIDE-VIS Sensitivity optimized for VIS at 500 nm	Qmini VIS/NIR
220 – 400 nm	370 – 750 nm	730 – 1080 nm	225 – 1000 nm	225 – 1000 nm	480 – 1100 nm
0.5 nm	0.8 nm	0.8 nm	1.5 nm	1.5 nm	1.5 nm